## Informace pro uživatele

CS

EN 420:2003+A1:2009; EN ISO 374-1:2016/Typ B; EN ISO 374-5:2016

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EN

EN 420:2003+A1:2009; EN ISO 374-1:2016/Type B; EN ISO 374-5:2016

Please read the instructions for use carefully

PPE of the III category

Complies with the requirements of Regulation (EU) 2016/425.

They are in compliance with the Regulation (EU) 2017/745 concerning medical devices (class I) and the Regulation of the European

Parliament and of the Council (EC) No. 1935/2004 on materials and articles intended to come into contact with food.

**Product:** Disposable Nitrile Gloves

<b>Design:</b>	STERN (4g) STERN ECO (3g) Non-sterile, Powder Free Nitrile Examination Gloves	<b>Code:</b>	3510 001 400 XX 3510 058 400 XX	<b>Size:</b>	S (6,5 - 7), M (7,5 - 8), L (8,5 - 9), XL (9,5 - 10)
----------------	--	--------------	------------------------------------	--------------	---

**Material:** NITRILE; no protein from natural latex, non-powdered / disposable, non-sterile.

**Purpose:** Handling dry and wet materials. The gloves provide their user with a time-limited protection against aggressive liquid chemicals – see the chemical table. Suitable for short-term contact with food, EN 1186. The gloves are intended for medical examination; they are a non-sterile medical item, class I, EN 455-1/2/3.

EN 420+A1 – protective gloves (general requirements)					Marking of gloves protecting against chemicals
(1-5)	Dexterity	Model class	5		
EN ISO 374-1 / Typ B					ISO 374-1/Typ B  KPT
Chemicals	Permeation Class (min. 0/max. 6)	Permeation time (minutes)	Degradation (%)	Appearance after chemical exposure	
K – 40 % Sodium hydroxide	6	> 480	- 36,4	No change	
P – 30% Hydrogen peroxide	2	> 30	18,5	Swollen	
T – 37 % Formaldehyde	5	> 240	14,1	Swollen	
EN ISO 374-5		LEVEL	Marking of gloves protecting against bacteria and fungi		
Protection against bacteria and fungi	pass		ISO 374-5:2016  VIBRI		
Protection against viruses	pass				

**Permeation time** for a given chemical product is the time after which the chemical permeates through the glove at the molecular level, sometimes without any visible damage to the glove. **Degradation index** is a degree of damage to the gloves when they are in contact with a specified chemical product; it manifests itself by change of their physical properties, e.g. softening, hardening etc.). The results are referred to as negative percentage degradation, when the glove material hardened after exposure to chemicals and required extra force when pierced.

**EN ISO 374-1:** This information does not reflect the actual duration of protection in the workplace and does not differentiate between mixtures and pure chemicals. The chemical resistance was only evaluated under laboratory conditions on samples taken from a palm (with the exception of 400 mm or longer gloves, when the cuff is also tested) and applies only to the tested chemicals. This resistance may vary if mixtures of chemicals are used. It is recommended to check that the gloves are suitable for the intended use, as workplace conditions may differ from the type test due to temperature, abrasion and degradation. In use, protective gloves may provide less resistance to hazardous chemicals due to changes in their physical properties. Moving, grinding, abrading, and degradation due to contact with chemicals, etc., can significantly reduce the actual usable time. Degradation can be the most important factor in choosing chemical resistant gloves for the purpose of working with aggressive chemicals. Prior to use, check the gloves for defects or imperfections (holes, cracks, tears...). If they are damaged, they must be disposed of. If gloves are damaged in any way during use, discard them immediately. **EN ISO 374-5:2016** Penetration resistance was assessed under laboratory conditions and applied only to the test sample.

**Attention:** Do not use damaged gloves. Do not use the gloves for handling sharp objects, flames or hot or heated objects. The gloves may only be used for the specified types and concentrations of chemicals. It is necessary to have the gloves tested in advance for each additional chemical. Protect the gloves against radiant heat and flames. It is the responsibility of every user to judge, whether the products are used in suitable conditions.

**Instructions for maintenance and decontamination:** After working with chemicals, do not touch the surface of the gloves with bare hands. Clean or wash the gloves properly before taking them off. Ensure so that no chemicals cannot penetrate through the cuffs. The gloves cannot be washed or chemically cleaned; they are intended for single use only.

**Storing and transport:** Gloves must be stored in a cool, dry place, out of direct sunlight. Avoid contact with petroleum products and grease. The lifespan of gloves depends on storage conditions. The shelf life is set at 5 years, provided the storage and transport conditions are observed. The month/year of production and expiration is marked on the box. Transport the product in the original packaging.

**Hazardous ingredients:** This product does not contain any substances known or suspected to adversely affect the hygiene or health of users. Nevertheless, the materials which may come into contact with skin of the user may cause allergic reaction at some sensitive people. Do not use gloves any more, if you notice any skin irritation or allergic reaction. In that case consult a doctor.

**Disposal:** Not used gloves shall be disposed of with household waste. After the gloves come into contact with chemicals, it is necessary to follow the relevant regulations for waste disposal of the respective chemicals. Treat them as chemical waste.

**Notified body** carrying out conformity assessment and periodic inspections: SATRA Technology Europe Ltd., Bracetown Business Park, Clonee, Dublin 15, D15 YN2P, Ireland. **NB: 2777.**

**Declaration of Conformity** see here: [www.canis.cz](http://www.canis.cz), in the bar of the individual products – “Documents for download”.

**Manufacturer:** PT Arista Latindo, Jl. Raya Bogor Km 36, Depok 16415 – Indonesia.

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MedNet GmbH Borkstrasse 10, 48163 Muenster, Germany.